

Help us test Ikigaro with your Ultrahuman ring

Ikigaro reads your blood panels and shows you your own trends over time. We have just built Ultrahuman syncing, and it has never once run against a real ring. Yours would be the first. About 10 minutes of setup, then you wear your ring exactly as you already do.

What we would read, and what we would not

Only the daily summaries Ultrahuman already computes: **sleep duration and score, HRV, resting heart rate, recovery, steps, blood oxygen, VO2 max and skin temperature deviation**. Plus, if you wear an M1 sensor, your daily glucose average, variability, time in target and estimated HbA1c.

We do not read your minute-by-minute glucose trace, your location, your workouts or anything else in your account.

Your data is shown only to you. It is never sold, never shared with anyone, and never used for advertising. Disconnect any time from Profile, which deletes our permission on the spot. We would only ever ask you to read numbers off your own screen and tell us whether they match.

What to do

1 Get in

Sign up at app.ikigaro.com with your email. You will land on a waiting list, which is normal: it is a closed beta. Tell us you have signed up and we will approve you within a few minutes.

2 Connect the ring

In the app, go to **Profile**, then **Connected devices**, and press **Connect** next to Ultrahuman.

Ultrahuman will ask you to approve access for *ikigaro*, listing Profile, Ring Data and CGM Data. Press **Approve**. You come straight back to Ikigaro.

3 Wear it for one full night

This is the part that cannot be rushed. Almost everything we want to check is measured while you sleep, so connecting and looking the same afternoon shows nothing.

4 The next morning, compare four numbers

Open your **Ultrahuman app** first and look at last night. Then open Ikigaro, go to **Profile > Connected devices**, press **Sync now**, then go to **Trends** and find the **From your devices** card.

Check these four against each other. They should match, and the point of the exercise is to find out if they do not.

In Ikigaro	Compare against, in the Ultrahuman app	What a mismatch would mean
Sleep , shown like 7h 05m	Your total sleep for the night	We may be reading their number in the wrong unit
HRV	Your overnight HRV, not the all-day figure	There are two HRV numbers and we may have picked the wrong one
Steps	Your step count for the same day	We may be reading an average instead of the daily total
Sleep score	Your sleep score for the night	Ultrahuman reports two different sleep scores

The failure we are actually hunting for

A number that is wrong but *looks* perfectly reasonable. If Ikigaro says 7h 05m and your ring says 7h 05m, that is a pass. If Ikigaro says 7h 05m and your ring says 6h 10m, that is exactly the bug we need, and nothing about the number would have looked suspicious on its own.

So please compare them side by side rather than judging whether each one seems plausible.

What to send back

A screenshot of the **From your devices** card in Ikigaro, and a screenshot of the same day in your Ultrahuman app. That is the whole report. If the four numbers match, say so and you are done.

Tell us too if something felt broken or confusing, even if it is not on this list. That is worth as much to us as the numbers.

If something does not work

No Ultrahuman option in Connected devices: your account is probably still on the waiting list. Tell us and we will approve it.

The card says nothing has synced, or shows no numbers: most likely the ring has not uploaded to Ultrahuman yet. Open the Ultrahuman app first, let it finish syncing your ring, then press Sync now in Ikigaro again. If it is still empty after that, tell us: that is a real finding and worth reporting.

No glucose numbers: expected, unless you wear an M1 sensor.

Thank you. This is the last untested piece of the integration, and one night of your ring's data settles it.

Questions, or to say you have signed up: reply to whoever sent you this, or write to team@ikigaro.com.