

Help us test Ikigaro with your ring or band

Ikigaro reads your blood panels and shows you your own trends over time. We have just built device syncing for **Oura**, **Ultrahuman** and **Whoop**, and it has never run against a real device. Yours would be among the first. About 10 minutes of setup, then you wear your device exactly as you already do.

What we would read, and what we would not

Only the daily summaries your device already computes: **sleep duration and score, HRV, resting heart rate, recovery or readiness, steps and blood oxygen**. On an Ultrahuman with an M1 sensor, also your daily glucose average, variability and time in target.

We do not read your workouts, your location, your minute-by-minute heart rate or glucose traces, or anything else in your account.

Your data is shown only to you. It is never sold, never shared with anyone, and never used for advertising. Disconnect any time from Profile, which deletes our permission on the spot. We would only ever ask you to read numbers off your own screen and tell us whether they match.

What to do

1 Get in

Sign up at **app.ikigaro.com** with your email. You will land on a waiting list, which is normal: it is a closed beta. Tell us you have signed up and we will approve you within a few minutes.

2 Connect your device

In the app, go to **Profile**, then **Connected devices**, and press **Connect** next to your device. Only devices we have finished setting up appear there, so if yours is missing, tell us.

Your device's own site will ask you to approve access for *Ikigaro* and list what it is granting. Press **Approve**. You come straight back to Ikigaro.

Whoop only: this needs an active Whoop membership, not just an account. If sign-in bounces you back without an approval screen, that is the reason, and it is nothing you did wrong.

3 Wear it for one full night

This is the part that cannot be rushed. Almost everything we want to check is measured while you sleep, so connecting and looking the same afternoon shows nothing.

4 The next morning, compare the numbers

Open **your device's own app** first and look at last night. Then open Ikigaro, go to **Profile > Connected devices**, press **Sync now**, then go to **Trends** and find the **From your devices** card.

Check these against each other. They should match, and the point of the exercise is to find out where they do not.

In Ikigaro	Compare against, in your device's app	What a mismatch would mean
Sleep , shown like 7h 05m	Your total sleep for the night, not your time in bed	We may be reading their number in the wrong unit, or counting time you were awake
HRV	Your overnight HRV, not an all-day figure	Most devices report two, and we may have picked the wrong one
Readiness	Oura: Readiness. Ultrahuman: Recovery. Whoop: Recovery	These are one number to us, so a mismatch means we read the wrong field
Steps	Your step count for the same day. Whoop does not count steps , so it should be missing there	We may be reading an average instead of a daily total

In Ikigaro	Compare against, in your device's app	What a mismatch would mean
Sleep score	Your sleep score. On Whoop, Sleep Performance	Some devices report two different sleep scores
Blood oxygen	Your overnight average SpO2, if your device measures it	This one was silently missing on Oura until recently

The failure we are actually hunting for

A number that is wrong but *looks* perfectly reasonable. If Ikigaro says 7h 05m and your device says 7h 05m, that is a pass. If Ikigaro says 7h 05m and your device says 6h 10m, that is exactly the bug we need, and nothing about the number would have looked suspicious on its own.

So please compare them side by side rather than judging whether each one seems plausible.

If you nap, that is a gift. Naps are logged separately from the night, and on two of these devices we found that an afternoon nap could quietly replace the whole night's sleep. A day where you napped is the single most useful day you can send us.

What to send back

A screenshot of the **From your devices** card in Ikigaro, and a screenshot of the same day in your device's own app. That is the whole report. If the numbers match, say so and you are done.

Tell us too if something felt broken or confusing, even if it is not on this list. That is worth as much to us as the numbers.

If something does not work

Your device is not listed in Connected devices: either your account is still on the waiting list, or we have not finished setting that device up. Either way, tell us.

Whoop bounces you back after sign-in, with no approval screen: Whoop only allows this for accounts with an active membership. Nothing you did wrong, and nothing we can fix from our side.

The card says nothing has synced, or shows no numbers: most likely your device has not uploaded to its own cloud yet. Open its app first, let it finish syncing, then press Sync now in Ikigaro again. If it is still empty after that, tell us: that is a real finding and worth reporting.

Some numbers are missing: often correct. Whoop does not count steps, not every device measures blood oxygen, and glucose only appears on an Ultrahuman with an M1 sensor. Tell us which are missing anyway, since that is exactly the kind of gap we are looking for.

Thank you. This is the last untested piece of the integration, and one night of your device's data settles it.

Questions, or to say you have signed up: reply to whoever sent you this, or write to team@ikigaro.com.